

## INQUIRY

When inquiry is made for parts or repair, please include the MANUFACTURING LOT NUMBER located on the handle. (See Figure 1).

## ACCESSORIES

CAT. #790 125cc HDPE Bottle w/3.65 tubing inlet and adapter  
CAT. #822 8-port disposable manifolds (25/box)

## REAGENT EXCLUSION LIST

This is not a complete listing — any reagent that is not compatible with Polypropylene or Viton® should be avoided.

## NOT TO BE USED

Amyl Chloride	Dipropylene Glycol	<b>SHORT TERM USE ONLY ONE DAY MAXIMUM, then rinse with water.</b>
Benzyl Alcohol	Ether	
Boric Acid	Ethyl Acetate	
Bromine	Hydrofluoric (HF) Acid	
Bromobenzene	Methyl Isopropyl Ketone	
Bromoform	Nitrobenzene	
Butadene	Oils: Cedarwood Cinnamon	
Butyric Acid	Oleum	
Carbon Disulfide	Perchloroethylene	
Chlorobenzene (Mono)	Tetrahydrofuran	
Chlorosulfonic Acid	Thionyl Chloride	
Chloroform	Toluene	
Chlorox	Trichloroethane	
Diethyl Benzene		
Diethyl Ether		

# StatMatic II

## 8-Channel Microplate Dispenser

### MODEL 980

## Operating Instructions

### SPECIFICATIONS

Dispensing Ranges - 0.8mL, 1.6mL, 2.0mL & 2.4mL (full-stroke)  
Dispensing per manifold port - 100uL, 200uL, 250uL, 300uL

Not Autoclavable

Only glass, Polypropylene, Platinum/Iridium, Teflon® and Ceramic contact the reagents (excluding flexible Tygon® Remote Inlet). See the recommended reagent exclusion list.

### INSPECTION

Carefully inspect your Dispenser for any shipping damage. Any defects or shortages should be reported to the supplier. Included with each Dispenser are the following accessories: (3) each of calibration rings, (1) bottle and cap assembly, (1) remote inlet tube assembly, (1) 8-port manifold.

### DISPENSE VOLUME SETTING (Refer to Figure 1)


All Dispensers are shipped set at the maximum dispense volume. Other dispense volumes may be set by installing the proper CALIBRATION RING. Unscrew the RETAINER CAP from the HANDLE. The PLUNGER BUTTON ASSEMBLY contains a RETURN SPRING and should be held firmly during unscrewing to prevent rapid separation and possible damage. Place the desired RING on the PLUNGER BUTTON and carefully guide the PLUNGER back into the opening in the GLASS BARREL.  
**(CAUTION:** Although the plungers are highly resistant to breakage, any severe bending or scratching may cause serious damage). Retighten the RETAINER CAP.

Teflon® Registered Trademark of Dupont



12555 Loma Rica Drive, Grass Valley, CA 95945 USA  
530-273-8888 800-WE PIPET (937-4738)  
FAX 530-273-2586

FIGURE 1



### SET UP

Included with the basic Dispenser is a 500mL Polyethylene Bottle with its Cap Assy. Connect the Cap Assy to one end of the long remote Inlet Tube. Connect the other end of this Tube to the side port of the Dispenser as shown in Figure 1. Install the 8-port Manifold to the end of the Dispenser as shown in Figure 1.

### PRIMING

Fill the Bottle with the desired reagent and remove the 8-port Manifold. Unscrew the Cap Assy from the Bottle and insert the outlet port of the Dispenser into the neck of the Bottle along-side of the tubing of the Cap Assy. Press and release the Plunger Button until bubble-free reagent is being discharged back into the Bottle. Reinstall the Cap Assy and the Manifold.

### OPERATION

For the highest order of precision in dispensing, the Plunger Button should be operated with consistent strokes. Avoid releasing the Button to “snap” back.

### CLEANING

The Dispenser should be cleaned at regular intervals, particularly if protein based reagents are used. One method calls for a suitable solvent or acid to be pumped through the Dispenser until it is purged of any residues. A second method requires the plunger to be removed from the barrel, which permits more thorough cleansing. (**CAUTION:** The plunger assembly contains a return spring and should be held firmly during unscrewing to prevent rapid separation and possible damage). Rotate the retainer cap counter-clockwise and carefully extract the Teflon® plunger. (Refer to Figure 1). (**CAUTION:** Although the Teflon® plunger is highly resistant to breakage, any severe bending will cause serious damage). Do not scrape or use abrasives to clean the plunger. A soft cloth and warm water are recommended. The barrel can best be cleaned with a small tube brush, followed by a thorough rinsing with water. Replace the return spring and carefully re-engage the plunger into the barrel. Rotate the retainer cap clockwise until fully seated.

If the Dispenser plunger appears to be stuck in the barrel, do not apply excessive force towards removal as plunger breakage may result. Place the entire unit in warm water and allow to soak for at least 10 minutes. Allow the Dispenser to cool. Then apply a gently twisting motion to the plunger. Generally this will free the plunger permitting removal for complete cleaning. Excessive removal pressure may cause a slight withdrawal of the barrel from its mounting flange. The proper, fully bottomed position is indicated by the handle spacer (Figure 1) being firmly in contact with the top of the barrel mounting flange. If displacement is noted, heat the flange in hot water and press downward on the upper portion of the handle until alignment is restored. Dispensers may be filled with distilled or deionized water for storage to minimize sticking of the plunger.

Teflon® Registered Trademark of Dupont