

# SPECIFICATIONS



## MiniPet

### VOLUMES

|      |      |      |      |       |       |      |
|------|------|------|------|-------|-------|------|
| 5 L  | 10 L | 11 L | 14 L | 20 L  | 25 L  | 30 L |
| 35 L | 40 L | 50 L | 75 L | 100 L | 200 L |      |

### PERFORMANCE (w/ generic Eppendorf style 0-200 L Tip)

| Volume                      | ±% Accuracy | ±% Precision (CV) |
|-----------------------------|-------------|-------------------|
| 100 L & 200 L (full volume) | 3           | 2                 |
| 75 L                        | 4           | 2                 |
| 50 L (full volume)          | 3           | 2                 |
| 35 - 45 L                   | 4           | 2                 |
| 20 - 30 L                   | 5           | 3                 |
| 10, 11 & 14 L               | 10          | 4                 |
| 5 L                         | 20          | 5                 |

### PHYSICAL

Dimensions (instrument only)  
3.75 in. Or 9.53 cm. Length  
.540 in. Or 1.37 cm. handle diameter

Weight  
.20oz. Or 6.25gm. (instrument only)

#### Appearance

|                |               |                    |                     |
|----------------|---------------|--------------------|---------------------|
| Red<br>5 L     | White<br>10 L | Pink<br>20 L       |                     |
| Orange<br>25 L | Plum<br>30 L  | Light Blue<br>35 L | Light Green<br>40 L |
| Yellow<br>50 L | Green<br>75 L | Blue<br>100 L      | Teal<br>200 L       |

### FLUID CONTACT MATERIALS

Tip Only (disposable)

### ACCESSORIES

0 - 200 L Yellow Tips  
DeTipper Bottle